

## Physics Experiment

### Absorption of $\beta$ -radiation in aluminum

WONG Man Hon, WONG Yuen Wah and ONG Chung Wo  
Department of Applied Physics, The Hong Kong Polytechnic University, Hung Hom, Kowloon,  
Hong Kong  
8 July, 2013

**A. Objective:** To study the absorption of  $\beta$ -radiation in aluminum (Al).

#### **B. Safety**

Never touch a radioactive source with bare hands. Whenever you need to transfer a radioactive material, always hold it with a lifting tools, like forceps. After the experiment, place the radioactive material back in a properly designed shielding case. Wash your hands thoroughly.

#### **C. Theory**

$\alpha$ -radiation can be stopped by a paper sheet, or air around 5 cm thick.  $\beta$ -radiation can be stopped by a Al sheet of around 5 mm thick.  $\gamma$ -radiation is much more difficult to be stopped (Fig. 1).  $I_0$  is the intensity of the incident radiation.  $I$  is the intensity of radiation after passing through an absorbing material with a thickness  $x$ .  $\mu$  is the linear attenuation coefficient (Fig. 2).

$$I(x) = I_0 e^{-\mu x} \quad (1)$$




Fig. 1 Penetration power of  $\alpha$ ,  $\beta$  and  $\gamma$ -radiation.




Fig. 2 Attenuation of radiation in a material of thickness  $x$ .

#### **D. Experiment**

##### **D.1 Safety assessment**

A survey meter show readings of \_\_\_\_\_ ( $\text{mR hr}^{-1}$ ) or \_\_\_\_\_ ( $\text{mSv hr}^{-1}$ ) for a source 30 cm from it. Compare the result with the data given in the document "Radiation units and safety". Hence, assess the risk.



### D.2 Radiation absorption

- Place a thin Al foil under the Geiger Muller (GM) tube (Fig. 3). Most of the  $\alpha$ -radiation is absorbed.  $\gamma$ -radiation is hardly detected by the GM tube. Therefore, the counts detected by the detector is mainly due to  $\beta$ -radiation.
- Turn on the setup by applying 400 volts to the GM tube.
- The number of counts from the background is collected for 300 s.
- Place a radium (Ra-) 226 isotope at the base plate of the setup. Increase the thickness of the Al sheets successively. Record the corresponding numbers of counts for 300 s for each thickness of Al. Subtract the background counts from the data and complete Table 1.
- Plot the net count rate against the thickness of Al.




Fig. 3 Setup for the measurements.

Table 1 Counts and count rates. Each measurement lasts for 300 s.

Thickness of absorber (mm)	0.00	0.212	0.371	0.602	1.21	2.12	3.14
Background counts (in 300 s)							
Background count rate (/s)							
Counts (in 300 s)							
Count rate (/s)							
Net count rate = Count rate – Background count rate (/s)							

### E. Analysis

- Determine the value of  $\mu$  for  $\beta$ -radiation absorption in Al from Eq.(1).
- According to the maximum annual effective exposure, discuss whether an Al container can be used to store a radioactive isotope emitting  $\beta$ -radiation.

- END -

